

# Experiment-3

## SUBCULTURE TECHNIQUES FOR PURIFICATION OF FUNGI

Primary fungal cultures often result in **mixed growth**. To isolate and purify a specific fungus, **subculturing techniques** are applied. These methods allow transfer of selected fungal parts (e.g., spores or hyphae) to fresh medium for **pure culture development**. The techniques used are:

- **A. Spotting**
  - **B. Plug Inoculation**
  - **C. Streaking**
- 

### Material

#### Common Materials for All Techniques

1. **Sabouraud Dextrose Agar (SDA) plates**
2. **Primary culture plate or old growth**
3. **Incubator (25±3°C for molds, 37°C for yeasts)**
4. **Parafilm or tape for sealing plates**
5. **Marker for labeling plates**
6. **Gloves, lab coat, and sterile workspace**

#### A. Spotting

7. Inoculating needle (platinum, aluminum, silver/iron types)
8. 70% ethanol for sterilization
9. Spirit lamp or Bunsen burner

#### B. Plug Inoculation

10. Sterile scalpel (optional)
11. Well borer (cork borer)
12. Forceps (optional)

#### C. Streaking

13. Inoculating loop
  14. Loop sterilizer or spirit lamp
-

## A. Spotting Technique


Used for filamentous fungi (molds)

Two variations:

- **Single Spot** – One central inoculation point
- **Multi-Spot** – Three points in triangular form

### Single Spot Technique – Procedure

1. Sterilize the inoculating needle:
  - If **platinum**: heat until red-hot.
  - If **aluminum** or **silver/iron**: dip in ethanol, then briefly flame.
2. Touch the sterile needle to a **colony from primary/old culture**.
3. Gently place it in the **center of a fresh SDA plate**.
4. Incubate upright at **22–28°C for 3–5 days**.



**Figure 3.1: Subculture by single spot technique**

## Multi-Spot Technique – Procedure

1. Sterilize needle as above.
2. Touch a colony from the primary culture.
3. Gently spot the spores onto **three points** on a fresh SDA plate in **triangle formation**, rotating the needle slightly at each point.
4. Incubate at **22–28°C for 3–5 days**.

💡 *Multi-spot technique provides better distribution of nutrients and reduces contamination risk compared to single spotting.*




Figure 3.2: Subculture by multi spot technique





Figure 3.3: Single spot of mold of left side and multispot of mold on right side

## B. Plug Inoculation

Used for transferring hyphal/mycelial plugs of molds

### Procedure

1. Sterilize the **well borer** by dipping in 70% ethanol, then flame to evaporate.
2. Press the sterile borer into a **primary culture plate** to extract a **mycelial plug**.
3. Transfer the plug to the **center of a fresh SDA plate**.
4. Gently press the plug to ensure contact with the medium.
5. Incubate at **25±3°C for 3–5 days**.



*Figure 3.4: Agar plug inoculation for mold*

---

## C. Streaking Technique


Used for yeast purification

Types of streaking:

- One-way
- Two-way
- Three-way
- **Four-way (recommended)**

### Procedure

1. Sterilize the **inoculating loop** and allow it to cool.
2. Touch a colony of yeast from primary/old culture.
3. Perform **four-way streaking** on an SDA plate for isolation.
4. Incubate **inverted at 37°C for 24–48 hours**.



**Figure 3.5: Streaking for sub-culture of yeast**



**Figure 3.6: Yeast growth on SDA**

## PRESERVATION AND REVIVAL OF FUNGI

To avoid repeating the lengthy process of isolation and purification, fungal cultures should be **preserved** in a way that maintains their viability and characteristics over time.

### 1. Sub-culturing


Short-term method; involves transferring fungi onto **PDA or SDA slants**.

#### Material

1. SDA/PDA slants
2. Inoculation tools
3. Refrigerator (4°C)

#### Procedure

- Inoculate slants and incubate:
  - **Molds:** 25°C for 7 days
  - **Yeasts:** 37°C for 48 hrs
- Store at 4°C
- **Revival:** Subculture onto fresh slants every 15 days



**Figure 16.1: Subculturing for revival of cultures**

## 2. Lyophilization (*For Yeast*)

### Procedure

- Place yeast culture in a **lyophilizer** to dry and powder
- **Revival:** Inoculate powder into broth and incubate under suitable conditions

*(Demonstration only)*

## 3. Cryopreservation

### Materials

1. Overnight broth culture
2. Sterile **glycerol (15%)**
3. Cryovials
4. -20°C or -80°C freezer

### Procedure

- Mix 1 mL broth culture with 150  $\mu$ L sterile glycerol
- Freeze at **-20°C or -80°C**

### Revival

- Thaw at room temp
  - Centrifuge at 6000 rpm (10 min)
  - Wash pellet with sterile saline, re-suspend, and inoculate on SDA
- 

#### 4. Soil Preservation (for molds)

##### Materials


1. Sterile soil in a bottle
2. Autoclave
3. Spore suspension (in normal saline)

##### Procedure

- Autoclave soil **twice** on consecutive days
- Add fungal spore suspension
- Incubate at **25°C for 10 days**
- Store at **room temp or 4°C**

##### Revival

- **Sprinkle soil** on SDA or
- Suspend in saline, settle, and spread supernatant on SDA



**Figure 16.2: Preservation of filamentous fungi in soil**

## 5. Microbank (Commercial Bead System)

### Materials

1. Cryovials with **polystyrene microbeads**
2. Buffer

### Procedure

- Transfer pure culture into vial
- Wait 2 min, discard buffer
- Microbes adhere to beads
- Store at **-20°C or -80°C**



## Revival

- Remove one bead
- Streak onto agar plate
- Incubate at recommended temperature

*(Demonstration only)*



**Figure 16.3: Cryovial with microbeads (left side) and microbank (Rightside)**

---